

## DETERMINATION OF PATULIN IN APPLE JUICE

**Regulations for apple juice:**  
 Europe (EC 1881/2006) : 50µg/Kg  
 USA (FDA CPG Sec.510.150) : 50µg/Kg  
**Regulations for apple juice for infants and young children:**  
 Europe (EC 1881/2006) : 10µg/Kg

### PROTOCOL OF PURIFICATION

#### Sample preparation

Loading solution: 2.5mL apple juice and 2.5mL of water-2% acetic acid are mixed.

#### Purification with a 3mL/100mg **AFFINIMIP® SPE Patulin** cartridge

##### Equilibration

- 2mL Acetonitrile
- 1mL water

##### Loading

- 4mL of loading solution

##### Washing of interferences (W1)

- 1mL NaHCO<sub>3</sub> in Water
- 2mL Water

##### Drying by applying vacuum 10 seconds

##### Washing of interferences (W2)

- 1mL Diethyl Ether

##### Elution (E)


- 2mL Ethyl Acetate

The elution fraction was then evaporated and dissolved in water containing 0.1% acetic acid before HPLC analysis.


#### HPLC Method

Column: Atlantis T3 column, 150mm x 2.1mm  
 Mobile phase: Deionized water/ACN (95/5, v/v)  
 Flow rate: 0.2mL/min  
 Detection: UV - 276nm  
 Injection volume: 100µL.

### RESULTS



Chromatograms of apple juice containing 25µg/kg of Patulin before (Red) and after (Blue) **AFFINIMIP® SPE Patulin** Clean-up



Chromatograms obtained after **AFFINIMIP® SPE Patulin** Clean-up of an apple juice spiked at 40µg/kg (tested twice, red) or at 10µg/kg (tested twice, blue) with Patulin or not spiked (orange)

Recovery of Patulin in apple juice after **AFFINIMIP® SPE Patulin** Clean-up and relative standard deviation calculated from results generated under reproducibility conditions.

Concentration of Patulin (ng/mL)	Recoveries %	% RSD <sub>R</sub>
10	97.9	11 (n=9)
40	90.6	11 (n=41)

#### Catalog number:

**3mL-100mg sorbent for most uses**  
 FS102-02 for 25 cartridges  
 FS102-03 for 50 cartridges  
**6mL-200mg sorbent for a DRIED APPLE or higher enrichment with apple juice**  
 FS102-02B -200mg for 25 cartridges  
 FS102-03B -200mg for 50 cartridges